


**DL 2101T60  
NURSE BOARD MODULE**

To be placed in the nurses room, it allows to stop the acoustic signal and to calm down the patient by sending a light signal to the patient room. It is composed of: a red calling lamp, a green accepted call lamp, a call buzzer, a light switch, a pushbutton and a bistable relay with two switching contacts.

Operating voltage: 24 Vac.




**DL 2101T61  
PATIENT ROOM BOARD MODULE**

To be placed in the patient room. It allows to call the nurse in case of need.

It is composed of: a red calling lamp, a green accepted call lamp, a red lamp outside the door, a pull button, a reset pushbutton and a bi-stable relay with 3 NO contacts.

Operating voltage: 24 Vac.



**DL 2100T  
WORKING BENCH WITH RACK**

Metal frame with three levels for fitting modules, assembled on bilaminated wooden table.

Square legs with compensator feet which can be adjusted to offset unevenness of the floor.

**DL 2100TA  
CABINET**

Fire varnished steel-plate, provided with shelves for module storage and key-locked doors.

Supplied with 4 rubber casters, the cabinet can be fitted under the DL 2100T bench.

